

# rubber dam



- Made of pure neutral rubber latex
- Vacuum-packed for prolonged shelflife
- Withstand autoclaving upto 118°C (224°F)
- Thickness: thin, medium or heavy
- Colours: green, dark or light
- In roll of 450 x 15,5 cm
- In cut pieces 5" x 5" 15,0 x 15,0 cm  
package of 52 pcs
- In cut pieces 6" x 6" 12,5 x 12,5 cm  
package of 36 pcs



— SWEDISH DENTAL SUPPLIES AB —

Södervägen 30, SE-232 52 Åkarp, Sweden  
E-mail: [info@swedent.se](mailto:info@swedent.se), Fax +46 40 46 05 13, Phone +46 40 46 02 01, [www.swedent.se](http://www.swedent.se)