


# screw posts

gold plated  
or titanium



Root canal screws conical shaped to avoid unwanted wedge effect. Electrolytically plated with a thick layer of 24k gold. Available in 6 sizes and 4 length. Screw posts in cylindrical or composite type are also available.

Size	1	2	3	4	5	6	
S	●	●	●	●	●	●	Short 8,0 mm
M	●	●	●	●	●	●	Medium 9,5 mm
L	●	●	●	●	●	●	Long 12,0 mm
XL				●		●	Extra long 14,0 mm
∅	1,05	1,20	1,35	1,50	1,65	1,80	Dimension

- ASSORTMENT BOX OF 60 PCS  
GOLD PLATED or TITANIUM  
3 pcs of 20 sizes, cross key and hallow key

- ASSORTMENT BOX OF 240 PCS  
GOLD PLATED or TITANIUM  
12 pcs of 20 sizes, cross key and hallow key

- ASSORTMENT BOX OF 120 PCS  
GOLD PLATED or TITANIUM  
6pcs of 20 sizes, cross key and hallow key

- REFILL BOX GOLD PLATED  
Box of 12 pcs one size

- REFILL BOX TITANIUM  
Box of 6 pcs one size



— SWEDISH DENTAL SUPPLIES AB —