




2000

temporary polycarbonate crowns




— SWEDISH DENTAL SUPPLIES AB —

Södervägen 30, SE-232 52 Åkarp, Sweden  
E-mail: [info@swedent.se](mailto:info@swedent.se), Fax +46 40 46 05 13, Phone +46 40 46 02 01, [www.swedent.se](http://www.swedent.se)

# TEMPORARY POLYCA




Art. no. 2195. Anteriors and Premolars in 60 different sizes. Ass. box of 180 crowns, with mould guides.



Art. no. 2190. Anteriors and Premolars in 60 different sizes. Ass. box of 180 crowns, without mould guides.

# CARBONATE CROWNS



Rib.



Art. nr. 2215. Refill box of 5 crowns of any mould or size.

## ORDERING INFORMATION

# MOULD GUIDES

Description	CROWN no						
Upper Right Centrals	14	13	12	11	10	101	100
Upper Left Centrals	15	16	17	18	19	102	103
Upper Right Laterals	24	23	22	21	20	2	
Upper Left Laterals	25	26	27	28	29	200	
Upper Right Cuspids	34	33	32	31	30	301	300
Upper Left Cuspids	35	36	37	38	39	302	303
Lower Anteriors (long)	64	63	62	61	60		
Lower Anteriors (short)	65	66	67	68	69		
Biscupids (large)	44	43	42	41	40		
Biscupids (small)	54	53	52	51	50		

The crowns are made of non-toxic polycarbonate. This plastic is due to its molecule structure extremely unsusceptible, even to the strain of important temperature changes and other physical or chemical influences. These very flexible and thin crowns can be cut, formed, bent, and pressed without breaking or ripping.

They form stable compounds with all acrylates and ZnO-tooth cements. Further these polycarbonate crowns are superior because of their good cervical fit, precise anatomy and shade stability.

### The following packages are available:

Art.nr.

- 2190     **Anteriors** and **Premolars** in 60 different sizes  
ass. box of 180 crowns **without** mould guides.
- 2195     **Anteriors** and **Premolars** in 60 different sizes  
ass. box of 180 crowns, **with** mould guides.
- 2215     **Refill box** of 5 crowns of any mould and size.
- 2220     **Loose crowns** in bulk, min 100 crowns each number.
- 2200     **Rib** with 7 crowns Upper Central Right or Left.
- 2205     **Rib** with 6 crowns Upper Laterals Right or Left.
- 2201     **Rib** with 7 crowns Upper Cuspids Right or Left.
- 2210     **Rib** with 5 lower Anteriors Long or Short.
- 2211     **Rib** with 5 crowns Biscupids Large (L) or Small (S)



## SWEDISH DENTAL SUPPLIES AB

Södervägen 30, SE-232 52 Åkarp, Sweden

E-mail: [info@swedent.se](mailto:info@swedent.se), Fax +46 40 46 05 13, Phone +46 40 46 02 01, [www.swedent.se](http://www.swedent.se)